

# PEDIATRIC(NEONATAL) OXYGEN CONCENTRATOR

MODEL: JAY-2



Pediatric  
(neonatal)  
machine



Small & light



0-2 L flow meter  
Minimum Scale 0.125L/Min



## Application

1. The flow meter of JAY-2 is 0-2L.
2. The minimum scale can be adjusted to 0.125L/min.
3. This model is designed specially for infants.

## Specification JAY-2

Flow Rate	1-2LPM
Purity	93% ( $\pm 3\%$ )
Outlet Pressure	0.03-0.07MPA
Sound Level	$\leq 45$ db
Power Consumption	$\leq 110$ w
LCD Display	Present working time, accumulated working time
Alarm	Power Failure alarm, High&low pressure alarm
Net Weight	9Kgs
Size	435×240×410mm

### Exclusive Distributor



**PT. PERMANA PUTRA MANDIRI**  
 MEDICAL EQUIPMENTS SUPPLIES AND SERVICES  
 Jl. Taman Mini Indonesia Indah (TMI)  
 Pintu II Atas No. 43 Cipayung, Jakarta Timur, Indonesia  
 Telp. : +62-21 841 5976, Fax: +62-21 840 7865  
 Email : info@ptppm.id, Website : www.ptppm.id