

BIN-4000A/4000B/4000C



Model	BIN-4000A	BIN-4000B	BIN-4000C
screen		LCD	
air mode		●	
skin mode		●	
humidity control		●	
O2 concentration control	-		●
scale	-	-	●
electronic lift		option	


LCD Touch Screen

- LCD touch screen displays temp./humidity/O₂ concentration
- 72 hours data record
- Various self-check alarms


Double Walls Hood


- Two side doors open
- Monitor tray


Integration Sensor Box

- Shows chamber temperature separately
- Independent over temperature alarm


Humidity Reservoir

- Removable for cleaning
- Humidity range:30%~90% RH


Baby Scale

- Under the bed, rang: 0.1~8KG

X-ray Tray under the bed
Specifications Infant Incubator

Model	BIN-4000A	BIN-4000B	BIN-4000C
Electrical requirement		~ 220V 50Hz	
Heater power		550VA	
Air mode temp.range		25.0°C ~ 37.0°C	
Skin mode temp. range		34.0°C ~ 37.0°C	
Air mode temp. alarm		< 38.0°C	
Skin mode temp. alarm		< 40.0°C (Air temp.)	
Air mode high temp. range		37.1°C ~ 39.0°C	
Skin mode high temp. range		37.1°C ~ 38.0°C	
High temp. mode alarm		< 40.0°C Air temp.	
Air mode temp. display range		20.0°C ~ 45.0°C	
Skin mode temp. display range		25.0°C ~ 45.0°C	
Air mode alarm		± 3.0°C	
Skin mode alarm		± 1.0°C	
Skin temp. sensor accuracy		± 0.3°C	
Temp. control accuracy		≤ 0.5°C	
Mattress temp. uniformity		≤ 0.8°C	
Humidity control range		30-90%RH	
Humidity display range		0-99.9%RH	
O ₂ cocentration range	-		20-60%O ₂
O ₂ cocentration alarm	-		± 5%O ₂
O ₂ cocentration display range	-		0-99.9%
O ₂ cocentration display accuracy	-		± 3%O ₂
Electrical scale weight range	-		0.1-8Kg
Electrical scale display accuracy	-		± 50g

Exclusive Distributor


PT. PERMANA PUTRA MANDIRI
 MEDICAL EQUIPMENTS SUPPLIES AND SERVICES
 Jl. Taman Mini Indonesia Indah (TMII)
 Pintu II Atas No. 43 Cipayung, Jakarta Timur, Indonesia
 Telp. : +62-21 841 5976, Fax : +62-21 840 7865
 Email : info@ptppm.id, Website : www.ptppm.id