




DISTRIBUTOR ALAT KEDOKTERAN INDONESIA
PT. PERMANA PUTRA MANDIRI
Jakarta, Indonesia

Telepon : 021 841 5976 / 021 877 93819
Email : marketing@alatkedokteran.id
website : www.alatkedokteran.id

ZONCARE'S PROGRESS DEPENDS ON YOUR TRUST

ZQ-660I

Full Digital Ultrasonic Diagnostic
System





Feature

- Original double wave beam imaging technology in the industry
- 10.1 inch medical use LED monitor
- 2 activated probe connectors, probe auto-freeze protection
- Efficient workflow and real-time operation navigation prompt makes the operation much more easier
- One key full screen display, hidden menu item, let the diagnosis more intuitive
- Measurement software package of Obstetrics, Urology, Gynecology, small Organ
- Data Port: USB2.0, Video, S-Video, RS-232, VGA
- Auto generate BMP format report, support USB import and export

Powerful Processing Analysis System



Report Generation



Complete Application



Full OB Formula



Various Body Mark



Patient Information



One Key Optimization



Color Images



3s Fast System boot




Full Screen View




Cutting Edge Digital Imaging Technology Realize Remarkable Ultrasonnd Image



Wide band frequency conversion technology



THI Imaging



High density and Broadband Frequency Transducers



Convex Transducer



Linear Transducer



Transvaginal Transducer



Micro convex Transducer