

Portable Ventilator MTV1000



- ✓ ICU Grade Ventilator
- ✓ Monitoring FiO_2 , options SpO_2 and EtCO_2
- ✓ Electronic PEEP and Blender control
- ✓ Built-in Air generator - Turbine blower type
- ✓ Built-in rechargeable Battery
- ✓ Dual & Single Limb circuit Compatible

MTV1000

Versatile Ventilator with Full Features

A variety of Ventilation Modes

MTV1000 is an ICU-grade ventilator designed for mobility. The essential functions and high performance of system can save your efforts for patient care.

Standard:

PACV, PSIMV, VACV, VSIMV, Spont+/CPAP, PRVC, Auto, Apnea back-up, O₂Stream®

Additional Options:

Pneumatic nebulizer (Requires nebulizer kit)

MTV1000 Highlights


- ✓ ICU, ER, transportable, general ward, medical care facility or silver hospital (Home care patient)
- ✓ ICU-grade performance ventilator with unique pneumatic technology
- ✓ Both Invasive and non-invasive ventilation applied.
- ✓ Dual and single limb circuit applicable

MTV1000 meets the essential requirements on performance, accuracy and high safety. The MTV1000 ventilator complies with:

- CE0470
- KFDA
- ISO13485
- ISO9001

Product Overview

UI



MTM1000 ICU Package



- 1 Alarm LED
- 2 7" Color Touch Screen
- 3 SpO₂, EtCO₂ Option modules
- 4 Exhalation Valve Assembly
- 5 Low Flow Oxygen Inlet
- 6 High Pressure Oxygen Inlet

MEKICS Ventilator



Distributor Alat kedokteran Indonesia
PT. Permana Putra Mandiri
Telepon : +62-021 841 5976
Fax : +62-021 840 7865
Email : info@ptppm.id
Website : www.ptppm.id

Jl. Taman Mini Indonesia Indah (TMII) Pintu II atas No. 43
Cipayung, Jakarta Timur – Jakarta, Indonesia

MEKICS Co., Ltd.
21, Sangjiseok-gil, Paju-si, Gyeonggi-do, 10911, Korea
T. 82 70 7119 2520 F. 82 70 5052 5800
sales@mek-ics.com
www.mek-ics.com