

# SH180 NCPAP



## Features

- 5.7 inch LCD display screen.
- Easy setting by common knob and touch key
- Advanced build in Air-Oxygen blender, easy setting
- Oxygen concentration from 21%-100%
- Independent flow meter
- Sensitive safety valves will prove safety to patients
- Will offer neonate intensive care
- Compatible with various generators



## Specifications

- ★ Ventilation Mode: CPAP
- ★ FiO<sub>2</sub>: 21-100%
- ★ Airway Pressure: 0-12cmH<sub>2</sub>O
- ★ Respiratory Rate: 0-60bpm
- ★ Flow Range: 0-15L/min
- ★ FiO<sub>2</sub> Upper Limit: 22-100%
- ★ FiO<sub>2</sub> Lower Limit: 21-99%
- ★ Airway Pressure Upper Limit: 0.1-12cmH<sub>2</sub>O
- ★ Airway Pressure Lower Limit: 0-11.9cmH<sub>2</sub>O
- ★ Alarm: Airway pressure high/low, FiO<sub>2</sub>



**PT. Permana Putra Mandiri**  
Medical Equipment Supplies and Services  
Jl. Taman Mini Indonesia Indah (TMII) Pintu II atas No. 43 Cipayung, Jakarta Timur - Indonesia  
Phone : 021 841 5976, Fax: 021 840 7865  
Email: [info@ptppm.id](mailto:info@ptppm.id), [marketing@alatkedokteran.id](mailto:marketing@alatkedokteran.id)  
[www.ptppm.id](http://www.ptppm.id) | [www.alatkedokteran.id](http://www.alatkedokteran.id)

