

300mA/125Kvp High Frequency Mobile X-Ray Unit

Model : WSM-300 / WSM-300P



300mA/125Kvp High Frequency Mobile X-Ray Unit

Model : WSM-300 / WSM-300P



SPECIFICATIONS

MODEL	WSM-300	WSM-300P
CONTROL UNIT	300mA at 100kVp	
HIGH VOLTAGE GENERATOR	200mA at 125kVp High Frequency Generator	
X-RAY TUBE UNIT	ROTATING ANODE TYPE	
TUBE SUPPORT UNIT	ELECTRO-MAGNETIC LOCK FOR MOVEMENT OF TRANSVERSE, VERTICAL AND TUBE ROTATION	
	LOCKING	
	VERTICAL TRAVEL	700mm ~ 1800mm
	TRANSVERSE TRAVEL	700mm ~ 850mm
	ROTATION OF VERTICAL COLUMN	±180° AT EACH 90° POSITIONS MANUALLY
ROTATION OF X-RAY TUBE UNIT	±180° AT EACH 90° POSITIONS	
ANGLE INDICATOR	CIRCULAR	
COLLIMATOR	MANUAL DRIVEN 24V / 100W MORE THAN 160Lux	
HIGH VOLTAGE CABLE	2.5m X 2 FEDERAL ENDS	
MOBILE CART	LOCKING SYSTEM WITH CASSETTE BOX	
POWER CONSUMPTION	SINGLE PHASE 220V, 50/60Hz 30KVA	SINGLE PHASE 220V, 50/60Hz 3KVA(PSU IN)
DIMENSIONS	MOBILE CART : 600W X 1,060L X 1,040Hmm / TUBE SUPPORT UNIT : 210W X 970L X 1,800Hmm	
NET WEIGHT	APPROX. 350KG	

Features

- Compact design & Color
- Single phase, High Frequency Inverter
- Easy operating procedure
- A.P.R Function internal organs



DISTRIBUTOR ALAT KEDOKTERAN INDONESIA
PT. PERMANA PUTRA MANDIRI
Jakarta, Indonesia

Telepon : 021 841 5976 / 021 877 93819
Email : marketing@alatkedokteran.id
website : www.alatkedokteran.id



ISO 13485:2003

DKQ-3504